

miScript™ miRNA PCR Array (384-well, 384HC)

Dog miRBase Profiler HC Plate 1 ptt

For mature miRNA expression profiling using SYBR® Green-based real-time PCR

Cyclers for use with miScript miRNA PCR Array formats

Format	Suitable real-time cyclers	Plate	Cat. no.
E	Applied Biosystems® models 7900HT (384-well block), ViiA™ 7 (384-well block); Bio-Rad CFX384™	384 well	331223 MIFD-3601AE
G	Roche® LightCycler® 480 (384-well block)	384 well	331223 MIFD-3601AG

Contents

miScript miRNA HC PCR Arrays Formats E and G consist of 2, 12 or 24 PCR array plates provided with 2, 12 or 24 Optical Adhesive Films.

Description

Profile the relative expression of up to 277 dog miRNAs using an exclusive PCR panel. The miScript Dog miRBase Profiler HC PCR array is a plate containing a total of 277 microRNA miScript Primer Assays. The set enables the accurate relative quantification of the most current dog miRNA genome (miRNome or microRNome) as annotated in miRBase (www.miRBase.org). This ensures scientists identify all of the miRNAs in body fluids, exosomes, cells, or tissues (fresh, frozen, or FFPE) that potentially impact a specific biological process. By testing the complete set of miRNAs, the miRBase Profiler enables the definition of superior miRNA biomarker panels. A set of controls also present on the array enables data analysis using the

For further details, consult the *miScript miRNA PCR Array Handbook*.

QIAGEN reserves the right to occasionally redesign individual assays on the miScript Arrays for improved performance. This revision history can be accessed by contacting technical support and supplying the batch numbers from your arrays.



Shipping and storage

miScript miRNA PCR Arrays are shipped at ambient temperature, on ice, or on dry ice depending on the destination and accompanying products. Upon receipt, store at -20°C . If stored under these conditions, miScript miRNA PCR Arrays are stable for 6 months after receipt.

Note: Ensure that you have the correct miScript miRNA PCR Array format for your real-time cycler (see table above).

Note: Open the package and store the products at the indicated temperature immediately on receipt.

Array layout

miScript Primer Assays can be ordered for each mature miRNA at www.qiagen.com/GeneGlobe.

Position	Mature miRNA ID	miScript Primer Assay Catalog #
----------	-----------------	---------------------------------

Related products

Product	Contents	Cat. no.
miScript II RT Kit (12)	For 12 cDNA synthesis reactions: miScript Reverse Transcriptase Mix, 10x miScript Nucleics Mix, 5x miScript HiSpec Buffer, 5x miScript HiFlex Buffer, RNase-Free Water	218160
miScript II RT Kit (50)	For 50 cDNA synthesis reactions: miScript Reverse Transcriptase Mix, 10x miScript Nucleics Mix, 5x miScript HiSpec Buffer, 5x miScript HiFlex Buffer, RNase-Free Water	218161
miScript SYBR Green PCR Kit (200)	For 200 reactions: QuantiTect SYBR Green PCR Master Mix, miScript Universal Primer	218073
miScript SYBR Green PCR Kit (1000)	For 1000 reactions: QuantiTect SYBR Green PCR Master Mix, miScript Universal Primer	218075
miRNome miScript miRNA PCR Array	Array of assays for the human, mouse, rat, or dog miRNome; available in 96-well, 384-well, or Rotor-Disc 100 format	Varies
miScript miRNA QC PCR Array	Array for assessing cDNA sample quality	Varies
miScript Primer Assay (100)	10x miScript Primer Assay (contains one miRNA-specific primer)	Varies
RT ² PCR Array Loading Reservoirs	12 x 5 ml capacity, irradiation-sterilized loading reservoirs	338162
384EZLoad Covers	Pack of 4 color-coded covers for loading 384-well plates	338125

For up-to-date licensing information and product-specific disclaimers, see the respective QIAGEN kit handbook or user manual. QIAGEN kit handbooks and user manuals are available at www.qiagen.com or can be requested from QIAGEN Technical Services or your local distributor.

Trademarks: QIAGEN®, miScript™ (QIAGEN Group); Roche®, LightCycler® (Roche Group); Applied Biosystems®, SYBR®, ViiA™ (Life Technologies Corporation); Bio-Rad®, CFX384™ (Bio-Rad Laboratories, Inc.).

1073798 08/2012 © 2012 QIAGEN, all rights reserved.

www.qiagen.com

Canada ■ 800-572-9613

Ireland ■ 1800 555 049

Norway ■ 800-18859

China ■ 8621-3865-3865

Italy ■ 800-787980

Singapore ■ 1800-742-4368

Denmark ■ 80-885945

Japan ■ 03-6890-7300

Spain ■ 91-630-7050

Australia ■ 1-800-243-800

Finland ■ 0800-914416

Korea (South) ■ 080-000-7145

Sweden ■ 020-790282

Austria ■ 0800/281010

France ■ 01-60-920-930

Luxembourg ■ 8002 2076

Switzerland ■ 055-254-22-11

Belgium ■ 0800-79612

Germany ■ 02103-29-12000

Mexico ■ 01-800-7742-436

UK ■ 01293-422-911

Brazil ■ 0800-557779

Hong Kong ■ 800 933 965

The Netherlands ■ 0800 0229592

USA ■ 800-426-8157



Sample & Assay Technologies